

HEPD-2

NANOPERSION PICCOLO

Ideal for laboratory use

- ❖ Improved dispersion performance by ultra-speed and stable operation.
- ❖ Can be used on a desktop since it is compact.
- ❖ Easy to operate.
- ❖ Removable mixing tank for easy cleaning.
- ❖ Power supply 100 volts.

Supports blending of any kind of material with high precision!



Pharmaceuticals



Feed



Cosmetics



Plastic



Metal



Ceramics





HEPD-2 High Efficiency Powder Dispenser NANOPERSION PICCOLO

Specifications

Model	HEPD-2
Tank capacity	1.5L
Jacket	Jacket type
Power supply	AC 100V single phase
Motor	0.6kW
Allowable speed range	16,000 rpm (Variable speed inverter)
Control Panel	Separate type with timer
Weight	20kg

External Dimensions



Kawata U.S.A., Inc.

712 W Algonquin Rd
Arlington Heights, IL 60005
Phone: 847-379-1449

www.kawata-usa.com

*These specifications are subject to change without notice.
*Please read the instruction manual carefully before use.